How to buy

- **1. Measure** your client to determine the size you need. (See dimension chart for help.)
- **2. Select** the appropriate Dynamic Stander and color.
- **3. Choose** the accessories you need.
- **4. Add** up the total cost.
- **5. Order** online in minutes at Rifton.com or see the last page of this catalog for other ways to order.
- **6. Call** us any time at **800.571.8198** if you need assistance. We're always here to help.



The **conversion kit** converts our smallest stander, the K110, to the next size up, the K130, for less than one quarter the cost of a new stander. Request from customer service.

Conversion kit K116 \$565

User dimensions (inches)	K110 mini	K130 small	K150 med.	K170 large
Height	25-32	30-40	37-52	48-65



Key user dimension: height

Select the appropriate Dynamic Stander by the user's overall height. Choose the model that allows for growth.

Important: User's weight must not exceed the maximum working load.

Item dimensions (inches)	K110 mini	K130 small	K150 med.	K170 large
Overall length and width	25 x 25½	25 x 25½	32 x 29	39 x 33½
Overall height	221/2-28	231/2-30	30-39	371/2-49
Width without big wheels	21	21	25	281/2
Wheel diameter	20	20	27	35
Width of body support	7–10	7–10	9-12	10-14
Height of body support	161/2-22	22-29	29-38	36-48
Item weight (lbs.)	28	30	40	65
Max. working load (lbs.)	50	50	120	175



